

Launching GenomeSuite Analyzer on EC2

To launch the GenomeSuite Analyzer AMI (Amazon Machine Image) using EC2, follow these steps:

Step 1: Start at the EC2 Dashboard

1. **Navigate to the EC2 Dashboard:**
 - Once on the EC2 Dashboard, confirm you are in the correct region.

Step 2: Select a Region

1. GenomeSuite Analyzer is available in four AWS regions:
 - * US East (N. Virginia) us-east-1
 - * US East (Ohio) us-east-2
 - * US West (Oregon) us-west-2
 - * Canada (Central) ca-central-1
2. Choose a Region:
 - * In the top-right corner of the EC2 Dashboard, click on the region selector drop-down.
 - * Select one of the four regions from the list of regions.

Step 3: Find the GenomeSuite Analyzer AMI

1. **Navigate to AMIs:**
 - In the left-hand menu under “**Images**”, click on “**AMIs**”.
2. **Filter for Public Images:**
 - At the top of the AMIs page, find the drop-down list labeled “**Owned by me**”. Click on it and choose “**Public images**”.
3. **Search for the GenomeSuite Analyzer AMI:**
 - In the search bar, type “**Sniffles_v.1.1.0**” and press Enter.
4. **Select the AMI:**
 - Find the AMI named “**Sniffles_v.1.1.0**” in the results.
 - Click the checkbox next to this AMI to highlight the row.
5. **Launch the Instance:**
 - With the AMI row selected, click on the “**Launch instance from image**” button at the top right of the page.

Step 4: Configure the Instance

1. **Name the Instance:**
 - On the “**Launch an Instance**” page, provide a name for your instance in the “**Name**” field.
2. **Select the Instance Type:**
 - Under “**Instance Type**”, select one of the following options:

For human exome, use these instance types:

1. p3.2xlarge
2. p3.8xlarge
3. p3.16xlarge

For human whole genome, use this instance type:

1. p4d.24xlarge

These instance types are optimized for high-performance computing tasks using Nvidia GPU's.

3. **Configure Key Pair:**

- Under “**Key Pair (login)**”, select the key pair you created earlier (e.g., "sniffles-keypair").

4. **Set Network Settings:**

- Under “**Network Settings**”, select “**Edit**”.
- Choose “**Select existing security group**” and select the security group you created earlier (e.g., "sniffles-sg").

5. **Configure Storage:**

- Under “**Configure Storage**”, set the storage size to **2048 GiB**.
- Ensure the storage type is set to **gp3** (General Purpose SSD).

Step 5: Launch the Instance

1. **Review the Settings:**

- Double-check all configurations to ensure everything is set correctly.

2. **Launch the Instance:**

- Click on the “**Launch Instance**” button at the bottom of the page.
- If you get an error message “Insufficient Capacity”, this means all of the available computing resources are in use.
At this point, you have two options:
 - 1) Wait a few minutes and try to launch the instance again.
 - 2) Choose a different Region and try to launch the instance again.

Typically, there are computing resources available in one of the regions.

3. **Instance Initialization:**

- Once launched, you will be taken to a page showing the status of your instance. The instance will take a few minutes to initialize.

4. **Verify the Instance:**

- After the instance is running, you can find it under “**Instances**” in the EC2 Dashboard. Make sure the instance status is “**running**” before attempting to connect.